

CoughAssist E70

Clearing airway secretions by simulating a real cough



Experiencing a
natural cough

Redefining non-invasive secretion clearance to improve patient lives

An effective cough is critical to keep airways clear

For those unable to cough or effectively clear secretions, deep suctioning is often used to clear the airway. Unfortunately, invasive suction methods can be uncomfortable for the patient and have been linked to complications such as hypoxia, tissue damage and infection.¹

The CoughAssist E70 provides an effective yet remarkably gentle non-invasive alternative for use in the hospital and at home.

CoughAssist therapy clears airways for longer periods of time than tracheal suctioning, and with fewer complications.²



Improving airway secretion removal

CoughAssist therapy has been clinically proven to increase peak cough expiratory flows and reduce recurrent respiratory infections.^{3,4}

The integrated Cough-Trak algorithm aids device **titration** and patient synchronization, helping both comfort and compliance. Adjustable oscillation levels enhance mobilization and increase the benefits of **therapy**.

CoughAssist E70 offers three customizable **therapy** setting presets to accommodate different patient conditions or circumstances once they have been **discharged** from the hospital.

Treatment integrated with the patient's life

An intuitive interface and large colour monitor make it easier to assess treatment and fine-tune device settings to improve **therapy** efficacy and comfort.

Device settings can be locked so that parameters cannot be inadvertently changed during treatment.

Mains or battery powered, the CoughAssist E70 is a truly portable solution, offering patients increased freedom and support.

Introducing new tools for close follow-up at home

Data management tools help assess therapy efficacy and adapt settings as required, or as a disease progresses.

- Peak Cough Flow and Tidal Volume are displayed after each cycle:
 - Tidal Volume monitoring helps to determine the proper inspiratory pressure needed to deliver a deep inhalation
 - Peak Cough Flow monitoring allows adjustment of the expiratory pressure needed to deliver an effective cough
- SpO₂ and heart rate monitoring at rest gives instant feedback on therapy efficacy
- An SD card records therapy data for extended follow-up
- Compatibility with DirectView software gives a complete view of therapy

CoughAssist E70 – for every step of the patient care pathway



Titration



Discharge




Home therapy



Follow-up

Cross platform continuity

Common design shortens product familiarization time while compatibility with our oximetry module allows monitoring of SpO₂. The detachable battery is also interchangeable with other Philips Respironics devices.



DirectView



Detachable battery and oximetry connectivity



Roll stand

Product information

Product information		
	Automatic mode	Manual mode
Preset	1, 2, 3	
Cough-Trak	OFF/ON	N/A
Inhale pressure	0 to 70 cm H ₂ O, in increments of 1 cm H ₂ O	
Inhale flow values	Low/Medium/High	
Inhale time	0 to 5 s, in increments of 0.1 s	N/A
Exhale pressure	0 to -70 cm H ₂ O, in increments of 1 cmH ₂ O	
Exhale time	0 to 5 s, in increments of 0.1 s	N/A
Pause time	0 to 5 s, in increments of 0.1 s Only if Cough-Trak is OFF	N/A
Oscillation	OFF/Inhale/Exhale/Both	
Frequency	1 to 20 Hz, in increments of 1 Hz. Only available if Oscillation is activated.	
Amplitude	1 to 10 cm H ₂ O, in increment of 1 cmH ₂ O. Only available if Oscillation is activated.	
AC voltage source	100 to 240 VAC, 50/60 Hz	
DC power source	12 VDC	
Dimensions (cm)	23.1 (h) x 29.2 (w) x 19 (d)	
Weight	3.8 kg (4.3 kg with battery)	

Ordering information

Device (part numbers)	International	Brazil	China	Japan
CoughAssist E70 Including device without battery, SD card, patient circuit large 1.8 m, carry bag, AC power cord, air filter and a circuit retainer	1098159	1098161	1098162	1098163

Accessories	Part number	Patient circuit kits	Part number	
			1.8 m tubing length	2.7 m tubing length
Oximetry interface cable	1098718	Patient circuit infant*	1090830	1098403
Foot pedal	1059017	Patient circuit toddler*	1090831	1098404
Roll stand	1098655	Patient circuit small*	1090832	1098405
Circuit retainer	1099035	Patient circuit medium*	1090833	1098407
Carry bag	1098884	Patient circuit large*	1090834	1098408
Water trap	1098720	Patient circuit trach**	1090835	1098409
Detachable battery	1043570	Patient circuit mouthpiece**	1090836	1098410

* contains: mask, tubing, mask adapter and bacterial filter

** contains: interface (mouthpiece or trach adapter), tubing and bacterial filter

- Invasive suction linked to complications such as hypoxia, tissue damage and infection.
AARC Clinical Practice Guideline. Endotracheal suctioning of mechanically ventilated adults and children with artificial airways.
Respiratory Care 1993;38(5):500-504.
- CoughAssist therapy keeping airways clear longer than trached suction and with fewer complications.
Sancho J, Servera E, Vergara P, Marin J. Mechanical in-exsufflation vs tracheal suctioning via tracheostomy tubes for patients with amyotrophic lateral sclerosis: a pilot study. Am J Phys Med Rehabil 2003;82(10):750-753.
- CoughAssist therapy clinical proven to increase Peak Cough Flows.
Chatwin M, Ross E, Hart N, Nickol AH, Polkey MI, Simmonds AK. Cough Augmentation with Mechanical Insufflation/Exsufflation in Patients with Neuromuscular Weakness. Eur Respir J: March 2003; 21(3):502-508.
- CoughAssist therapy proven to reduce recurrent respiratory infections.
Alice C. Tzeng and John R. Bach. Prevention of Pulmonary Morbidity or Patients with Neuromuscular Disease. Chest 2000;118: 1390-1396. DOI 10. 1378/chest 118.5.1390.

The CoughAssist E70, delivering **innovations** to meet your patient's needs



Lightweight, portable and yet robust

Detachable battery delivering
1 day of therapy*.
Handle for easy transportation.



Large colour screen and intuitive interface

Customizable according to
environment (detailed/limited
view).

Displays instant feedback
(Peak Cough Flow, Tidal Volume
and SpO₂).

Flexibility in delivery of therapy

Control lever for manual
application of therapy.

Foot pedal so that the caregiver
can do manual chest thrust while
holding the interface.

Automatic mode, with
Cough-Trak trigger option.



Data Management

SD card records more than one
year of therapy data that can be
queried via DirectView software.

Enhanced therapy efficacy and comfort

Cough-Trak feature for initiation
of therapy by the patient
triggering inspiration.

Oscillation to help enhance
mobilization.

3 presets of settings available.

* 1 day of therapy is defined as performing a typical treatment 4 times.
A typical treatment being 4 to 6 sequences of 4 to 6 cough cycles at +/-40 cm H₂O.

For more information please contact us at:

Air Liquide Healthcare, Stratus House, Unit 1 IDA College Park, Blanchardstown Road North, Blanchardstown, Dublin 15, D15 PEC4.

Tel: +353 (0)18091800 Freephone: 1800 24 02 02 (ROI only) 0800 328 5875 (NI only)

Fax: +353 (0)1829 3966 Email: healthie@airliquide.ie Website: www.airliquidehealthcare.ie