

Respireo Primo P

simple & ergonomic pillow mask



We are all different

The Respireo Primo series is a range of masks for ventilation that has been specifically designed to fit and adapt to the individual needs of each patient.

Simple & ergonomic, Respireo Primo P is a pillow mask combining silence and ease of use in order to quickly win the compliance of patients.

Available in vented and non vented versions, it is the reference mask for all patients willing to keep their liberty and autonomy in their daily activities.

It is also recommended for patients inconvenienced by forehead support, contact on the nasal bridge or suffering from claustrophobia.

Respireo
Primo P

Respireo Primo P Mask



Adaptable tubing

Ball joint system for a rotation at 360°.


Adjustable fixation: over the head, on the lateral sides or free downwards.

Clear visual field.



Stable headgear

Designed to ensure comfort and stability. Two levels of adjustment.



Ergonomic cushion

Airtight and fitting the nostrils thanks to the large sides and the 'bellows' technology.

Easy assembly of the cushion with a visual and tactile poka-yoke.



Silent and diffuse leak

Ordering information

Product information	Size	Part number
Respireo Primo P		
Respireo Primo P (Manufacturer reference 472105*)		KM219900
Respireo Primo P Non Vented (Manufacturer reference 472106*)		KM220000
Respireo Primo P Cushion		
Respireo Primo P Cushion (Manufacturer reference 472210)	Small	KM220700
Respireo Primo P Cushion (Manufacturer reference 472211)	Medium	KM220800
Respireo Primo P Cushion (Manufacturer reference 472212)	Large	KM220900

Respireo Primo P is also available in a **non-vented version** for patients ventilated with active valve ventilators.



(*) Medical device class IIa 0051.

For more information please contact us at:

Air Liquide Healthcare, Stratus House, Unit 1 IDA College Park, Blanchardstown Road North, Blanchardstown, Dublin 15, D15 PEC4.

Tel: +353 (0)1 809 1800 Freephone: 1800 24 02 02 (ROI only) 0800 328 5875 (NI only)

Fax: +353 (0)1 829 3966 Email: healthie@airliquide.ie Website: www.airliquidehealthcare.ie