

AirMed TRAVEL-AIR

portable compact compressor



AirMed TRAVEL-AIR

AirMed TRAVEL-AIR is a low cost option for a portable compact compressor; designed to efficiently deliver prescribed medication to the bronchial lung passages. Developed for the successful treatment of asthma, allergies and other respiratory conditions.

Key features

- Portable hand-held compressor
- Can be used worldwide
- Rechargeable battery for mobile use included
- Comes complete with a starter pack, AC adapter and travel bag
- Compact and lightweight

AirMed TRAVEL-AIR

portable compact compressor

Dimensions	106 x 66 x 128mm (w/o battery) 106 x 66 x 166mm (with battery)
Weight	395g (w/o battery) 712g (with battery)
Electrical requirements	12VDC from AC adapter or battery pack
Flow range and pressure range	5.0 - 6.7 L/min, 8 - 10 psi 50 - 70 kPa, 0.5 - 0.7 bar
Particle size	MMAD 4 microns
Avg nebulisation rate	>0.20ml/min
Compliance	Electrical Safety Class II Type BF. Continuous
Battery operating time	45 minutes
Accessories included	Neb cup x 1, driveline x 1, adult face mask x 1, child face mask x 1, mouthpiece x 1, connector x 1, inlet filters x 5, universal AC adaptor x 1, rechargeable battery x 1, travel bag x 1, rubber base x 1

Specifications are subject to change without notice.

Product	Part Codes
AirMed TRAVEL-AIR (UK plug)	3605902
AirMed TRAVEL-AIR (EU plug)	3605903
Adult Year Pack	3605545
Child Year Pack	3605544
Car Adaptor	3605918
Spare Rechargeable Battery	3605917



For more information please contact us at:

Air Liquide Healthcare, Stratus House, Unit 1 IDA College Park, Blanchardstown Road North, Blanchardstown, Dublin 15, D15 PEC4.

Tel: +353 (0)1 809 1800 Freephone: 1800 24 02 02 (ROI only) 0800 328 5875 (NI only)

Fax: +353 (0)1 829 3966 Email: healthie@airliquide.ie Website: www.airliquidehealthcare.ie