

F&P AIRVO™ 2

humidifier with integrated flow generator



AIRVO™ 2

Humidifier with integrated flow generator.

The AIRVO™ 2 generates high flows of warmed and humidified respiratory gases, delivered to the patient through a variety of nasal, tracheostomy and mask interfaces.



AIRVO™ 2 Features

With its integrated flow generator and innovative oxygen delivery system, the AIRVO™ 2 provides an effective, versatile solution for the delivery of oxygen therapy.

High Performance

- Designed to exceed the ISO standard for respiratory humidification
- Up to 60 L/min of supplemental oxygen


Versatile

- Variety of Optiflow patient interfaces for upper and bypassed airways
- Wide flow range: 2-60 L/min

Set up and use

- Graphical display menu for settings
- Integrated O₂ mixing
- Inbuilt O₂ sensor
- No probes or external air supply required

Accessories



AIRVO™ 2 specifications

Product description

AIRVO™ 2 Humidifier with Integrated Flow Generator

Intended use	The AIRVO™ 2 humidifier is for the treatment of spontaneously breathing patients who would benefit from receiving high flow warmed and humidified respiratory gases. This includes patients who have had upper airways bypassed. The flow may be from 2 - 60 L/min depending on the patient interface. The AIRVO™ 2 is for patients in hospitals and long-term care facilities
Humidity output	Complies with ISO 8185:2007 standard
Flow settings	10 to 60 L/min (in increments of 5 L/min) (default) 2 to 25 L/min (in increments of 1 L/min) (Junior mode)
Patient interfaces	F&P Optiflow™ nasal cannulas F&P Optiflow™ tracheostomy interfaces Masks with standard 22 mm medical taper
Models	PT101AZ / PT101EE / PT101EW / PT101UK / PT101US
Dimensions	295 mm x 170 mm x 175 mm (11.6" x 6.7" x 6.9")
Weight	2.2 kg (4.8 lb) - unit only 3.4 kg (7.5 lb) - packaged in bag incl. accessories
Electrical ratings	50-60 Hz 100-115 V ~ 2.2 A (2.4 A max) 220-240 V ~ 1.8 A (2.0 A max)

For more information please contact us at:

Air Liquide Healthcare, Stratus House, Unit 1 IDA College Park, Blanchardstown Road North, Blanchardstown, Dublin 15, D15 PEC4.

Tel: +353 (0)1 809 1800 Freephone: 1800 24 02 02 (ROI only) 0800 328 5875 (NI only)

Fax: +353 (0)1 829 3966 Email: healthie@airliquide.ie Website: www.airliquidehealthcare.ie