

VisionAire™ 5

Stationary Oxygen Concentrator



Time-proven design
Patient satisfaction
Ease of service


CAIRE®
O₂ is what we do.



The foundation of your patients' oxygen therapy

Compact and lightweight oxygen delivery

The VisionAire™ 5 has been completely redesigned to be state-of-the-art, meeting all current compliance standards with quality engineering and reliability.

Time-proven design

- **Industry-leading reliability**—greater than 99%* reliability rate out of the box
- **Comprehensive warranty**—3-year warranty on the device
- **Proven durability**—engineered with robust valves and sieve beds to minimize failure rates

Promotes patient satisfaction

- **Effective oxygen delivery**—up to 95.5% concentration from 1.0 to 5.0 LPM, meeting a wide variety of patient prescriptions
- **Quiet**—you will hardly know this lightweight and compact device is there
- **Energy efficient**—reduced power consumption to minimize the burden on the patient's energy bill

Ease of service

- **Maintenance free**—no required preventative maintenance reducing your total cost of ownership
- **Dual hour meter**—clear LCD provides time of use in hour and minutes
- **Enhanced diagnostics**—on board alarm alerts facilitate the ability to quickly troubleshoot

Technical specification

Flow Rate	1.0 – 5.0 lpm
Weight	13.6 kg
Size (l x w x h)	528 x 358 x 292 mm
Noise level	39 dBA
Power	290 W (all flow rates)
Maximum Tubing	15.2 m (50 ft)
Warranty	3 years



*Data on file.

Serving your patients across the continuum of CAIRE.

Non-Delivery Solution

Reduce operating expenses and time in the field while improving patient care with the CAIRE nondelivery solution—the combination of the industry-leading FreeStyle® Comfort® portable oxygen concentrator with the maintenancefree VisionAire™ 5 stationary oxygen concentrator.



O₂ is what we do.

For more than 50 years, CAIRE has been designing and delivering the best technologies and services to support oxygen therapy. CAIRE products are trusted globally by patients, providers, clinicians and the military. Let us put our experience to work for you.

For more information please contact us at:

Air Liquide Healthcare, Stratus House, Unit 1 IDA College Park, Blanchardstown Road North, Blanchardstown, Dublin 15, D15 PEC4.

Tel: +353 (0)18091800 Freephone: 1800 24 02 02 (ROI only) 0800 328 5875 (NI only)

Fax: +353 (0)1829 3966 Email: healthie@airliquide.ie Website: www.airliquidehealthcare.ie