

DreamStation Go

Travel PAP

Anywhere you go



DreamStation Go convenience without compromise

Whether you're traveling for business or venturing out on the vacation of a lifetime, DreamStation Go provides an easy, reliable and portable PAP experience.

- Our smallest and lightest tubing¹
- Our smallest PAP machine
- Our lightest PAP at less than two pounds
- Complies with FAA requirements²

Packs performance anywhere you go

Wherever your travels take you, sleep apnea follows. That's why DreamStation Go is an indispensable travel partner, packs the same performance as our home CPAP machines into a light, streamlined package. It even has a convenient USB port to charge your mobile device.

Learn more

Contact your provider today to learn more about DreamStation Go.

www.philips.com/dreamstationgo

¹ When compared to Philips 22mm performance tubing

² RTCA/DO-160G section 21, category M.



Choose the mask you want



Mask compatibility

DreamStation Go gives you the freedom to choose the mask you want!¹ There's no need for you to switch from the mask you are currently wearing at home when using a travel CPAP.

CPAP therapy even 'off the grid'



Overnight battery

DreamStation Go's optional overnight battery can help you stay committed to your therapy when you are off the grid. DreamStation Go, set at 10cm of pressure displayed an average run time of 13 hours² when running with the DreamStation Go overnight battery.

Simple therapy adjustments



Intuitive colour touch screen

Intuitive colour touch screen navigation doesn't require a mobile device app to adjust your therapy.

¹ Masks with standard 22mm outer diameter connection.

² Battery run times may vary due to higher CPAP pressure settings, excessive mask leakage, altitude, temperature, charging devices with the USB port, and other variables.

Innovation for people on the go

DreamStation Go provides the portability and packability travelers want and the effective therapy they need.



Charge

Only travel PAP with integrated USB port allows for easy charging of mobile devices.



Power

Built-in power supply reduces the number of parts for you to manage for added convenience.



Safe

Complies with FAA requirements for in-flight use.¹



Connected

Bluetooth is standard with every device.



Control

No app or mobile device needed to adjust your therapy settings.



Performance

Device is compatible with our DreamMapper monitoring self-management tool so you can continue to track your progress.



¹ RTCA/DO-160G section 21, category M.

For more information please contact us at:

Air Liquide Healthcare, Stratus House, Unit 1 IDA College Park, Blanchardstown Road North, Blanchardstown, Dublin 15, D15 PEC4.

Tel: +353 (0)18091800 Freephone: 1800 24 02 02 (ROI only) 0800 328 5875 (NI only)

Fax: +353 (0)1829 3966 Email: healthie@airliquide.ie Website: www.airliquidehealthcare.ie