


Wisp

rediscover dreams



Paediatric nasal mask



Minimal contact, maximum comfort

The comfort and performance of the Wisp minimal contact nasal mask is specifically designed for children. Wisp paediatric applies a child-friendly fabric pattern, a modified cushion curvature, and family support tools to provide a positive experience for your smallest patients.



Unique design for unique patient needs

The cushion, frame and headgear are designed to fit the sizing and bone structure of paediatric patients. The open line of sight allows children to comfortably wear the mask during their normal daily activities, while the Leak Correction Dial allows parents to easily make adjustments without removing the mask.



Simplified ordering for increased efficiency

Wisp paediatric features a streamlined product offering to help make it easier to stock and supply to your customers. Give your entire paediatric population the comfort they need with Wisp Paediatric and Wisp youth.

Paediatric Nasal Mask

Cushion design

Featuring auto seal groove technology, delivering minimal contact and a snug seal.

Click fit headgear clips

With adjustable straps for quick and easy adjustments.

360° rotating elbow

Gives the child freedom to move.

Leak Correction Dial

For less intrusive mask adjustment. If leaks occur, turn the dial left and right, gently nudging the cushion closer to the face.

Optional tube wrap

A tubing management loop enables the tube to be positioned over the head and away from the front of the body.

Product information

Part number

Wisp Paediatric nasal mask

Fit pack with 3 cushion sizes (SC Small, SC Medium, SC Large) with child size headgear

1104961

Wisp Youth nasal mask*

Fit pack with 3 cushion sizes (Petite, Small/Medium, Large) with adult reduced size headgear

1109308

* Intended for both adult and Paediatric use (7 years old or >18 kgs).



For more information please contact us at:

Air Liquide Healthcare, Stratus House, Unit 1 IDA College Park, Blanchardstown Road North, Blanchardstown, Dublin 15, D15 PEC4.

Tel: +353 (0)1 809 1800 Freephone: 1800 24 02 02 (ROI only) 0800 328 5875 (NI only)

Fax: +353 (0)1 829 3966 Email: healthie@airliquide.ie Website: www.airliquidehealthcare.ie