

Patient's name:	Address:	
Date of birth:		
Diagnosis:	Eircode:	
Next of kin:	Email address:	
NOK telephone number:	Private patient? Yes <input type="checkbox"/> No <input type="checkbox"/>	
Patient's telephone number(s):	Medical card number:	Expiry Date:
Weight: *weight in KG is required for neonatal, infant and paediatric patients only		

Device selection	
I wish to order a spot check finger clip device <input type="checkbox"/> (go to signatures section)	I wish to order a continuous monitoring device <input type="checkbox"/> (go to continuous monitoring section)

Continuous monitoring (not needed for spot check device)	
Please select probes Adult finger clip: <input type="checkbox"/> Paediatric finger clip: <input type="checkbox"/> Reusable sensor: <input type="checkbox"/> Disposable sensor: <input type="checkbox"/> <small>(recommended for neonates, infants and paediatrics)</small>	If using a reusable probe Frequency of bandage changes: <input type="text"/> <small>(supplied in a box of 50/100)</small> Frequency of probe changes: <input type="text"/> <small>(supplied in a box of 24)</small> One probe will be provided free of charge at installation
Alarm settings Patient's oxygen saturation alarm limits: Upper: <input type="text"/> Lower: <input type="text"/>	Patient's heart rate alarm limits: Upper: <input type="text"/> Lower: <input type="text"/>

Additional details

Signatures (this form must be signed by hand)	
Print name:	Hospital:
Signed:	Ward & room number:
Position:	Bleep number:
Contact telephone number & email (should we need to contact you to clarify settings):	

Send to
Medical card patients: Fax to HSE PCCC for approval. Requests to a PCCC office must include a current medical card number. Non medical card patients: Fax direct to Air Liquide Healthcare on 01 829 3966.

Accessory information

Accessory	Weight limitation	Qty	Notes
Paediatric finger clip	>5kg <45kgs	1	
Adult finger clip	>15kgs	1	
Reusable neonatal probe	<3kgs	10	Only plasters need to be reordered
Reusable infant probe	3-15kgs	10	Only plasters need to be reordered
Disposable neonatal probe	<3kgs	24	Only recommended for paediatrics
Disposable infant probe	3-20kgs	24	Only recommended for paediatrics
Disposable paediatric probe	10-50kgs	24	Only recommended for paediatrics
Disposable adult probe	>30kgs	24	

Nonin 7500 specifications

- For continuous use
- Convenient, lightweight & compact
- Fits easily on a tabletop or into its own carrying case
- Versatile: Monitors neonate to adult patients in acute care to homecare environments
- Easy to use: Bright LED displays and large display screen offer clear readability
- Durable: Rugged design
- Powerful: Over 70 hours of memory available for extended monitoring
- Oxygen Saturation Display Range: 0-100% SpO₂
- Pulse Rate Range: 18-321 BPM



Apex spot check device specifications

- Can be used with adults and children
- Measures SpO₂ 70-99%
- Pulse rate 30-235 BPM
- Turns off after 8 seconds of being idle to conserve battery
- 2xAAA batteries provide up to 30 hours use



If you require any further information before you complete this form, please contact our customer service team on 1800 24 02 02, and ask to speak with one of our clinical nurse specialists.