



**NEW**



For treatment of diseases  
of the lower airways



**PARI COMPACT2** Inhalation Device  
Quality at an affordable Price

## Nebuliser Therapy for the Lower Airways:

### **PARI COMPACT2** Inhalation System

PARI stands for fast, effective treatment. As a result of more than 110 years experience as specialists in nebuliser therapy, PARI delivers devices of the highest quality. Quality at an affordable price.

- Robust & reliable compressor device
- Easy to clean
- Includes the clinically-proven PARI LC PLUS nebuliser
- The clinical performance and safety of the PARI LC PLUS nebuliser used for different inhalation solutions was proven in more than 90 clinical studies.
- Recommended nebuliser in the SmPCs of 15 inhalation solutions for diseases of the lower airways such as COPD, asthma & bronchitis.
- Made in Germany – 100% quality control

## The advantages of nebuliser therapy

- No coordination necessary
- Can be used with all medications approved for inhalation
- Varying doses possible depending on medication used
- Additional moisturisation of the airways to promote mucociliary clearance



## PARI COMPACT2 Inhalation System

- PARI LC PLUS nebuliser
- PARI child mask soft



	Item No.		IFU Languages
PARI COMPACT2 Inhalation System	152G1001	Western Europe	de, en, es, fr, it, nl
PARI COMPACT2 Inhalation System	152G1006	Eastern Europe	cs, el, et, hr, hu, pl, ro, sk, sl

Useful Accessories	Item No.
PARI LC PLUS Nebuliser	022G8104
PARI PEP S System	018G4000
PARI LC Interrupter	022G1000
PARI Filter/Valve Set	041G0500

Aerosol Characteristics:	
Total Output Rate:	440 mg/min
MMD:	3.6 $\mu$ m
Mass fraction < 5 $\mu$ m:	67 %
<i>Measurement with Helos SympaTec at 23 °C and 50% relative humidity. Nebulised medium: 0.9% NaCl (5 ml). Inspiratory flow 20 l/min.</i>	

# PARI SERVICE CENTER

+49 (0)8151 279 220

## We're here to help!

PARI's International Service Coordinators in Starnberg, Germany, offer an expert service and are dedicated to assisting you personally with any questions regarding:

- Use and handling of PARI inhalation devices
- Cleaning and servicing
- Accessories and spare parts
- Product availability and lead-times
- Order processing and logistics

If you need information about a repair or service for a PARI device, please visit **[www.pari.com/partners](http://www.pari.com/partners)** to find a distributor and/or service partner in your country.

We look forward to hearing from you.

Simply give us a call or send an email to **[export@pari.de](mailto:export@pari.de)** for a prompt response.



PARI GmbH  
Moosstraße 3  
82319 Starnberg  
GERMANY  
Email: [export@pari.de](mailto:export@pari.de)  
Tel. +49 (0)8151 279 0  
Fax +49 (0)8151 279 101  
[www.pari.com](http://www.pari.com)