

# WellNeb

## Portable Compressor Nebuliser

WellNeb is a compressor nebuliser engineered for ultimate portability, perfect for those on the go. Its compact size, lightweight design, and convenient carrying handle ensure it remains accessible at all times.

The fine particle sizes and reliable output rate provide efficient treatment, quickly alleviating symptoms for users. WellNeb's user-friendly products offer simple operation, assisting in managing health, and maintaining a high quality of life.



# WellNeb

## Portable Compressor Nebuliser

- WellNeb is a compressor nebuliser designed for easy portability, especially for people on the move.
- Easy filter replacement, anti-drool mouthpiece.
- The compact size, lightweight and convenient carry handle of the WellNeb facilitate portability and allow users to keep the WellNeb accessible at all times.
- Motor-piston drive ensures reliable, consistent performance with quiet operation.



### User Friendly Design

The notched design provides secure placement for the air tub.



### Effective and Efficient Treatment

Delivers aerosolised particles MMAD: 5µm, Dispenses 4ml of medication in 20 minutes.



### Lightweight and easy to carry

The convenient handle allows for easy portability and mobility, with a lightweight design that weighs only 1.36kg.



### Accessories

According to your region.

## Specifications

Item	Specification
Power	220V, 50 Hz
Power Consumption	180W
Medication Capacity	6mL
Particle Size	0.5µm to 10µm
MMAD	5 microns
Sound Level	Less than 60dB at 1 metre
Average Nebulization Range	≥0.2ml / min
Compressor Pressure Rating	> 25 psi
Operating Pressure Range	8 - 13 PSI
Litre Flow Range	5.0 - 8.0 Lpm
Operating Temperature Range	5°C to 40°C (41°F to 104°F)
Operating Humidity Range	15% RH - 90% RH
Storage Temperature Range	-25°C to 70°C (-13° to 158°F)
Storage Humidity Range	0% RH - 90% RH
Dimension (L x W x H)	238 mm x 178 mm x 99 mm (9.37" x 7.01" x 3.90")
Weight	1.36 kgs
Standard Accessories	Nebuliser, air tube, mouthpiece, filters. Optional Accessories: Adult mask, child mask

Specifications are subject to change for improvement without notice

Wellell Inc.  
Print-2024.10 All rights reserved

